

Jefferson Radiology Image Exchange Instructions

Click the email link from a shared study and you will be brought to the below page. If you have created an account in the past, click 'Already have an account?' and log in. If you have not registered, fill out the required information and read and accept the terms.



Register

[Already have an account?](#)

For instructions and information on using Jefferson Radiology's Study Share please [click here](#)

* First name

* Last name

* Email Address

* Password

* Confirm password

8 characters or longer, containing 2 numbers and 2 letters

* I accept the [Terms of Use](#)


* I accept the [Privacy Policy](#)

* I accept the [Indications for Use](#)

Register

Cancel

Once you are logged into your account, your screen will look like the below with a list of the studies that have been shared with you.



The screenshot displays the Jefferson Radiology web interface. At the top left is the logo, and navigation links for 'Studies' and 'Messages' are visible. On the top right, there are links for 'Help', 'Lisa Castle', and 'Sign Out'. The main section is titled 'Studies' and includes a '+ Add a Study' button. A table lists a study with the following details:

Study ID	Study Name	Date/Time	Actions	Share
TEST ^ TEST (F) DOB 04-02-1986 MRN X107634ECHN	BONE ^ LOWER EXTREMITY PR 6 images Acc#: RPhys: Sarker ^ Gopal	08-31-2021 12:29 PM	Image Reports Download	Share

Below the table, there is a pagination control showing '100 items per page' and navigation arrows. At the bottom right, it indicates 'showing 1 study detail report'.

On this page you can view the images in the web viewer by clicking *images* or upload them into your system by clicking *download*.